

## Release Note Version 3.3.3

UPDATE  
3.3



SKILLQUBE®  
Training

# qubeSERIE Update 3.3.3

Released on 10.11.2021

## qubeSERIE

### Software is ready for iOS 15

Apple recently released iOS 15 for all their devices. We checked the compatibility of our software with the new operating system and can now recommend you to update your devices to the newest version

### Don't update immediately

With the release of any new operating systems by Apple we don't recommend updating your devices immediately. Always wait for our official approval so we can ensure that our software still runs perfectly. We always notify you if you are good to update your devices!

### Optimierungen:

- Optimized Translations
- Bugfixes and improved performance

# UPDATE 3.3



**SKILLQUBE®**  
Training

## **NEW RELEASE**

### **NEW Media center integrated in all qubeSERIES monitors**

-Receive media now directly on the monitor, e.g. 12-lead ecg.

### **IMPROVEMENT qubeCPR**

-Now it will be easier to identify the qubeCPR sensors, with the option to assign colours and numbers via the qubeCONTROLLER.

### **IMPROVEMENT:**

- Optimizations to the sound engine
- Optically modified login dialog
- Fixed a bug where SpO2 could increase to 101%
- Optimization of the SpO2 waveform for a realistic display of pressure fluctuations
- Optimizations for the support area
- General performance improvements

qube<sup>ZERO</sup>



### NEW qube<sup>ZERO</sup>:

The latest member of the qubeSERIES is now available. This versatile generic simulation monitor comes with a new refreshing screen design and new functions

- Switch between different views including clinical, numerical and prehospital
- 12 lead ECG
- Pacer
- AED and manual defibrillator
- CPR feedback integration
- Metronome in resuscitation mode
- Media center

---

### qube<sup>CONTROLLER</sup>

NEW:

- Portrait mode has been deactivated
- The app cannot be exited accidentally
- Sensor calibration can now be canceled
- In the preview, Media can now also be zoomed into

Bug fix:

- Pressing analysis deactivates the artifacts and thus resolves a UX issue



---

### Get in touch with us:

Mail: [support@skillqube.com](mailto:support@skillqube.com)  
Service-Hotline: +49 (0)6222 – 38 62 65  
Web: [www.skillqube.com/support](http://www.skillqube.com/support)

## qube<sup>3</sup> / qube<sup>3T</sup>

### NEW:

- Charging tone for the qube3 / qube3T was added
- Function extension for the chart while in CPR mode
- Automatic detection of lead II or pads

### Optimization:

- EtCo2 is now also displayed in DEFI / CPR mode
- Voice instructions in qube3 / qube3T optimized

### Bug fix:

- Correction of errors with PADS and synchronized cardioversion



## qube<sup>15</sup>

### Optimization:

- If SpO2 is activated, the heart rate is no longer displayed

### Bug fix:

- Incorrect functionality fixed (loading the defibrillator using the rotary push wheel)



## qube<sup>7</sup>

### NEW:

- Adjustments to SOFT8 from SCHILLER
- PACER adapted to SOFT8

### Optimization:

- The alarm LED now also lights up when there is an alarm
- Adjusted display in manual defi
- Blood pressure measurement display adapted



## qube<sup>X</sup>

### NEW:

- Integration of the IABD measurement

### Optimization:

- Revised functionality when switching on
- Adjustments when opening the pacer menu



## qube<sup>Assessment</sup>

### Optimization:

- Reference values are now also displayed in laboratories
- Optimizations for iOS 12.5

