



SKILLQUBE

Quick Start Guide

qube**SERIES**



skillqube.com




Controller iPad

Monitor iPad

Assessment iPad

1 Setting up devices






1. Establish an Internet-capable WiFi connection on all iPads.
2. Next, open the **Settings** app  > **Sign in to your iPad** to sign in to your devices with your Apple ID. If you don't have an Apple ID, you can set one up in the same menu.
3. After logging in with your Apple ID on all devices, open the **App Store** app and search for **SKILLQUBE**.
4. Download the **qubeCONTROL** app for the controller and the **qubeASSESS** app for the assessment. Download the interface available to you for the monitor, e.g., the **qubeZERO** app.



2 Activate devices

1. With the delivery of the new simulation device, the access data for our **qubeCLOUD** were sent to the device administrator. Please follow the instructions in the e-mail to create your user account. You will need that access data in the next step to activate the license.
2. Now open **qubeCONTROL**. When you open the app for the first time, the software requires you to log in to our **qubeCLOUD**. Enter the device administrator access data from step 1.
3. **qubeCONTROL** will now search for available licenses. Select the available license to activate your device. Activation may take a short time because the first scenarios and device-specific settings are already being transferred.
4. In the next step, open **qubeASSESS** to log in with your **qubeCLOUD** access data. Tap on **Not signed in** in the top right corner.

By logging in, existing scenario-based media will be loaded into the assessment, reducing the amount of data during the simulation. This step is not mandatory, but we recommend it.



3 Connect devices



Before you start the simulation, ensure your devices are on the same WiFi network. Then open **qubeCONTROL**, **qubeASSESS**, and the desired monitor.

1. To connect controller and monitor, tap the red flashing button on the top right corner (of the controller): **Connect Device**.
2. Tap to open a view of all available nearby devices. Locate your monitor and tap to select it.



Need help with using the iPad?
For more information, visit
support.apple.com


3. You will be asked to confirm the connection on the monitor. If the devices are successfully connected, a blue frame will appear.
4. Repeat this connection procedure with the assessment and, if available, the ventilation simulator.
5. To start a scenario, tap the top red flashing bar **Choose scenario**. You can now choose between a free scenario or a pre-programmed scenario. Once you have selected the desired scenario, tap **Load** at the bottom of the screen and then **Start now**.



4 Running a simulation

After starting the scenario, you can activate the vital signs and adjust various settings.


1. Activate and deactivate vital signs by simply tapping on them. Depending on the vital signs, you need to make further settings. If the button for the vital sign value, e.g., ECG, is red, the value is active and visible on the monitor.
2. The active vital signs can be changed using the sliders. To modify the waveform, tap the desired waveform at the bottom.



3. To send the changed vital signs to the monitor, tap **Send** in the lower right corner of the controller. You can use the +/- buttons above to set a trend time over which the vital sign value should develop.
4. To end a pre-programmed scenario, click the **pause** symbol at the top right. To end a scenario, tap the top red flashing bar **Pause or end scenario**. You can view the logbook immediately after the scenario. The simulation monitor switches off automatically.



Note

Internet connection: An Internet connection is required to activate the license and to install updates for the **SKILLQUBE** apps and the Apple operating system (iPadOS). We recommend using the included WiFi router to ensure a stable device connection during the simulation training. The access data can be found on the router. Connect all iPads to this network. To disconnect known WiFi networks, go to the **Settings** app  **WiFi > i > Forget This Network**.

General information

Keep your software up to date! Subscribe to our newsletter to receive information about software updates for our apps.

In addition, other simulation products are available to extend the **qubeSERIES**, such as our **qubeAIRFLOW**, **qubeCPR** and **qubeTRIGGER**.

All information can be found on our website:
skillqube.com



Notes

Contact person

Activation date

Other

Contact

SKILLQUBE GmbH

Hohenaspen 20
69231 Rauenberg
Germany

+49 6222 38 62 65
support@skillqube.com
skillqube.com/support