

## Release Note Version 3.6



UPDATE

3.6

qubeSERIE

# qubeSERIE Update 3.6

Release Date October 4th, 2022

### General improvements:

- The device will check if local networks are active. Local networks are required to connect the iPads
- Performance Improves and bugfixes

## qubeCloud

- Error descriptions optimized
- Complemented ventilation and intensive care values
- App settings got moved to qubeCLOUD
- Define pager texts in qubeCLOUD
- Create Lab reports in qubeCLOUD
- PDF Bibliothek, match files to specific scenarios
- Layout designer for qubeC3 and qubeC3T



---

## qubeCONTROL

- Reengineering for Ventilator Control Interface
- Improved synchronization for scenarios and media files
- Synchronize paper texts and Lab Reports with qubeCLOUD
- Print ECGs from the Media Center on your controller
- New ECGs available
- Improved Phase Description for a better view



## qube3 & qube3T

- During AED analysis Metronom will be deactivated
- Create your monitor layout by utilizing the qubeCLOUD Layout designer



---

## qubeAED

- Optimized voice prompt
- Implemented a variety of standard AED scenarios for stand-alone training



---

## qubeMS2

- New ventilation modes: PCV, BiLvl ASB, CPR IPPV, Man. CPR
- New views available featuring flow and pressure waveforms



# qubeASSESSMENT

- Assessment will now show pre-implemented SOPs from qubeCLOUD



## Contact us for questions

Mail: [support@skillqube.com](mailto:support@skillqube.com)

Service-Hotline: +49 (0)6222 – 38 62 65

Web: [www.skillqube.com/support](http://www.skillqube.com/support)