

Update 3.0.0 qubeAVgo


Release Notes Version 3.0.0 qubeAVgo

Note: iOS update required!

To use the latest version of **qubeAVgo**, you must install at least iOS 17.

Highlights in this update

- Setup assistant
- Selectable storage location



qubeAVgo

- Setup assistant supports you when setting up a new video recording
- USB image sources can now be integrated (*USB camera not possible at the same time as iPad camera)
- Integrate HDMI video sources using the HDMI grabber
- Integrate the microphone via the iPad's jack port
- You can now change the storage location of videos; offline and iCloud storage are possible
- Camera control improvements
- Video export optimizations
- Language adjustments and optimizations
- Optimized and restructured user interface for a better overview
- Improvements in the timeline view
- Logout in the app is now possible (*no effect on license binding on iPad)

Contact:

support@skillqube.com

+49 6222 – 386265

skillqube.com/support