

Update

6.2.0

qube**SERIES**

Release Notes Version 6.2.0

General

- Stability and performance optimizations
 - **qubeBOX** update to version 1.2.3
 - Improved display of connection quality
 - Optimizations in ventilation simulation
-

New interface


- **qube3000** — simulates the Dräger Oxylog 3000 plus

qube**3000**



qubeCONTROL

- Realistic simulation with a dedicated selection of electrodes and sensors



- Completely modernized CPR dashboard with improved analysis and feedback functions



- Updated logbook
- New notification system for providing relevant status and system information

SKILLQUBE

- Connection indicator (Wi-Fi or peer-to-peer): The connection indicator shows whether the devices are connected directly to each other via Wi-Fi or peer-to-peer and whether the connection is active and stable.



qubeC3 & qubeC3T

- Advanced options in the layout editor

qubeAED

- Metronome support during the CPR phase

Important note

- When upgrading to version 6.2.0, please ensure that you perform the update on all iPads in use.
- The change also requires an update of the **qubeBOX** firmware to version 1.2.3. To do this, all sensors must be disconnected during the update. Instructions for the update can be found [here](#).

Contact
support@skillqube.com
+49 6222 – 386265
skillqube.com/support